




CQ-03140-D2O-DEPT135-271009-298 K-500 MHz

DEPT135-mmcd-clem D2O /opt/topspin ravi {1 A4 - 104}



NAME CQ-03140-D2O-271009  
EXPNO 5  
PROCNO 1  
Date\_ 20091029  
Time 6.34  
INSTRUM spect  
PROBHD 5 mm CPTXO 13C  
PULPROG dept135  
TD 65536  
SOLVENT D2O  
NS 256  
DS 16  
SWH 29761.904 Hz  
FIDRES 0.454131 Hz  
AQ 1.1010548 sec  
RG 203  
DW 16.800 usec  
DE 30.00 usec  
TE 298.0 K  
CNST2 145.0000000  
D1 2.00000000 sec  
D2 0.00344828 sec  
D12 0.00002000 sec  
TD0 1

===== CHANNEL f1 =====  
NUC1 13C  
P1 10.00 usec  
P2 20.00 usec  
PL1 3.00 dB  
PL1W 36.25705338 W  
SFO1 125.6974365 MHz

===== CHANNEL f2 =====  
CPDPRG2 waltz16  
NUC2 1H  
P3 15.70 usec  
P4 31.40 usec  
PCPD2 100.00 usec  
PL2 1.25 dB  
PL12 17.33 dB  
PL2W 19.96668053 W  
PL12W 0.49238619 W  
SFO2 499.8419994 MHz  
SI 32768  
SF 125.6848680 MHz  
WDW EM  
SSB 0  
LB 1.00 Hz  
GB 0  
PC 1.40